

> Innovative products  
that help people.

**reha medi**  
Hoffmann GmbH

# Mobilizer Norbert



Discover more



**ERGONOMIC, SAFE  
AND COMFORTABLE**



 developed. designed. built.  
made in Germany



> Product video on:  
[www.fruehmobilisation.com](http://www.fruehmobilisation.com)

[www.rehamedi.de](http://www.rehamedi.de)



## Product features

The patient remains in a stable position when the backrest is tilted, shear forces on the patient are reduced and the nursing staff are relieved by minimising positioning and repositioning measures.

An optional height adjustment, mechanical or electrical depending on the selected equipment, ensures the setting of an ergonomic working height for the nursing staff and provides the positive psychological aspect: the patient is at eye level with his environment.


The viscoelastic and removable upholstery with wipe-disinfectable covers relieves pressure and enables longer, gentle sitting, even for immobile patients.

A wide range of colours is available. Swivelling armrests allow barrier-free transfer even with the back in an inclined position. This is particularly helpful for patients with respiratory insufficiency.

Suitable additional equipment is available for almost every application. We will be happy to advise you on your choice.

**> We will be happy to advise you!**

## Equipment options Examples



Art.-Nr.	6040106	6040116	7040115	6040110	6040109	6040107
Designation	Headrest fixing	Positioning pad	Leg strap	Chest strap	Lap belt	Calf strap

### > Technical data

Reclining height adjustable ..... 58 - 80 cm  
 Forward tilt ..... 80°  
 Standing position ..... 90°  
 Head-down position ..... 10°  
 Length of the lying surface ..... 190 cm  
 Seat width ..... 57 cm  
 Total width ..... 74 cm  
 smooth running wheels ..... Ø 150 mm  
 own weight ..... ca. 90 kg  
 Maximum permissible patient weight ..... 150 kg

#### > Widened version

Seat width ..... 65 cm  
 Total width ..... 79,5 cm

#### > Extension option up to 250 kg

Maximum permissible patient weight ..... 250 kg

### > Colour selection

CAPRI lemon	CAPRI lime	CAPRI kiwi
CAPRI orange	CAPRI rot	CAPRI lavendel
CAPRI azur	CAPRI delft	CAPRI anthrazit

