

> Innovative products
that help people.

 **reha medi**
Hoffmann GmbH

Curalizer



Mehr entdecken



easy to clean



Power+
optional



**IDEAL FOR CARE, TRANSPORT
AND PATIENT MOBILISATION**



 developed. designed. built.
made in Germany



Product video on:
www.fruehmobilisation.com

www.rehamedi.de

Product features

The patient remains in a stable position when the backrest is tilted, shearing forces on the patient are avoided and the nursing staff are relieved by minimising positioning and repositioning measures. A position-independent height adjustment, mechanical or electrical depending on the selected equipment, ensures the setting of an ergonomic working height for the nursing staff. The viscoelastic and removable upholstery with wipe-disinfectable covers is pressure-relieving and enables

longer, gentle sitting. A wide range of colours is available. Swivelling armrests allow barrier-free transfer even when the back is tilted - particularly helpful for patients with respiratory insufficiency. The removable head section extends the possibilities in patient care. Four swivelling, large castors with central brake system and directional lock ensure safe and easy manoeuvrability even in confined spaces.

Equipment options



Art.-Nr. 8040300
Designation neck cushion



8040100
Table with cushion



5040005
harness



6040110
Chest strap



6040110
Sitting trousers



8040701
Storage roll



Art.-Nr. 6040119
Designation Stabilisation cushion



6040202
Infusion stand



6050001
Additional battery

>

Factory options



8240300
Widening

> Technical data

Seat height adjustable.....	45 - 73 cm
Seat depth.....	45 cm
Seat width.....	50 cm
Total width.....	60 cm
Overall length (sitting position).....	95 cm
smooth running wheels.....	Ø 125 mm
own weight.....	ca. 65 kg
Maximum permissible patient weight.....	230 kg

> Widened version Art.-Nr. 8240300

Seat width.....	60 cm
Total width.....	71 cm

> Colour selection

